



## Zimmer® One-Piece Implant



### An esthetic solution from a trusted source.

Introducing the *Zimmer* One-Piece Implant, an exciting new addition to the Zimmer Dental family. Engineered to perform as well as the trusted and successful *Tapered Screw-Vent*® Implant, *Zimmer* One-Piece Implants bring greater simplicity and convenience to the restoration process and deliver the esthetics your patients desire. The simple one-piece design allows for a single-stage protocol, saving time and resulting in less trauma to soft tissue, and less discomfort for your patients.

#### Zimmer One-Piece Implants feature:

- One-piece construction, eliminating micromovement and maximizing strength
- An integral, pre-contoured abutment, saving time
- Zimmer's proprietary *MTX*™ surface, allowing for increased bone apposition and promoting osteoconductive capacity<sup>1,2</sup>
- A tapered body in a variety of sizes, enabling placement in all areas

<sup>1</sup> Todisco M, Trisi P. Histomorphometric evaluation of 6 dental implant surfaces after early loading in humans. *J Oral Implantol* 2006;32. In press.

<sup>2</sup> Trisi P, Marcato C, Todisco M. Bone-to-implant apposition with machined and *MTX* microtextured implant surfaces in human sinus grafts. *Int J Periodontics Restorative Dent* 2003;23(5):427-437.




## No compromises.

Zimmer One-Piece Implants feature contoured abutments with prepared margins, so they deliver an integral abutment that functions like a prepared tooth, with little or no additional adjustment required. Available in a range of diameters and lengths, and in both straight and angled versions, a Zimmer One-Piece Implant is available for virtually all tooth locations. And since the design is based on the trusted *Tapered Screw-Vent* Implant, they're a smart, convenient solution that can be immediately provisionalized and/or immediately loaded.\*

- The tapered implant body is similar to a natural tooth root, expanding placement options while allowing for initial implant stability
- Titanium alloy construction and multi-lead threads offer maximum strength and stability
- 3.0mmD Zimmer One-Piece Implants allow placement in the limited tooth spacing of the maxillary laterals and mandibular lateral and central incisors
- 3.7mmD and 4.7mmD Zimmer One-Piece Implants are indicated for placement in all tooth locations
- Snap-on impression caps and a complete assortment of other contoured copings simplify the restorative process

\* Not for 3.0mm Zimmer One-Piece Implants.



### Zimmer One-Piece Implants, Straight, with MTX Surface

Implant Diameter	Abutment Emergence Profile	10mmL	11.5mmL	13mmL	16mmL
3.0mmD	3.5mm	ZOP30S10	ZOP30S11	ZOP30S13	•
3.7mmD	4.5mm	ZOP37S10	ZOP37S11	ZOP37S13	ZOP37S16
4.7mmD	5.5mm	ZOP47S10	ZOP47S11	ZOP47S13	ZOP47S16

### Zimmer One-Piece Implants, 17°, with MTX Surface

Implant Diameter	Abutment Emergence Profile	10mmL	11.5mmL	13mmL	16mmL
3.7mmD	4.5mm	ZOP37A10	ZOP37A11	ZOP37A13	ZOP37A16
4.7mmD	5.5mm	ZOP47A10	ZOP47A11	ZOP47A13	ZOP47A16

### Zimmer One-Piece Surgical Kit, Complete – ZOPTST

To learn more about Zimmer One-Piece Implants, please visit us online at [www.zimmerdental.com](http://www.zimmerdental.com), or to speak to a sales representative, call 800 854 7019.

To receive our eNews visit us at [http://www.zimmerdental.com/news\\_eNewsletterSignUp.aspx](http://www.zimmerdental.com/news_eNewsletterSignUp.aspx)

For more information about our Products, Professional Programs and Educational Opportunities, contact us:

In the U.S. 800 854 7019  
 To fax an order 888 225 2483  
 Outside the U.S. +1 760 929 4300  
 Australia +61 (0)2 9950 5444 or 1 800 241 916  
 Canada +1 905 567 2073 or 1 800 265 0968

France +33 (0)1 45 12 35 35  
 Germany +49 (0)761 4584 722/723  
 Israel +972 (0)3 6124242  
 Spain +34 93 846 05 43



Zimmer Dental  
 1900 Aston Avenue  
 Carlsbad, CA 92008-7308  
 USA

[www.zimmerdental.com](http://www.zimmerdental.com)